



Achilles

VACOped[®] by OPED



Indications

- Achilles tendon rupture
- Ankle fractures
- Metatarsal fractures
- Calcaneus fractures
- Severe distortions
- Arthrodesis

The VACOped and VACOped short are cast-stable function orthosis for the immobilization of the foot and ankle.

With the ability to provide a range of motion, the VACOped allows for greater freedom of movement than a plaster cast and has been proven to shorten a patient's recovery process.

- The most comfortable and effective boot option for Achilles ruptures
- Option of range of motion (ROM) from -15° to +30° for VACOped and -15° to +20° for VACOped Short
- Vacuum cushion forms perfectly to individual patient 's anatomy
- Changeable liners for hygiene and comfort
- Easy wound inspection simplifies treatment
- Removable sole for better hygiene
- Moveable toe plate to lengthen the plantar footbed

OPED
medical inc.

5545 St. Jacques O. Montreal, QC H4A 2E3
O: 514.369.3311 F: 514.369.1911
Toll Free: 1.877.222.3311 Fax Toll Free 1.877.527.1911
www.bmiortho.com

BMI
BERRETTA MEDICAL INC.